

# HSPB6 Antibody, HRP conjugated



PACO51619

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

cytoplasm, nucleus, protein homodimerization activity.

**Gene ID:**

HSPB6

**Uniprot**

O14558

**Synonyms:**

Heat shock protein beta-6 (HspB6) (Heat shock 20 kDa-like protein p20), HSPB6

**Immunogen:**

Recombinant Human Heat shock protein beta-6 protein (1-160AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A