

HSPB6 Antibody, Biotin conjugated



PACO51621

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytoplasm, nucleus, protein homodimerization activity.

Gene ID:

HSPB6

Uniprot

O14558

Synonyms:

Heat shock protein beta-6 (HspB6) (Heat shock 20 kDa-like protein p20), HSPB6

Immunogen:

Recombinant Human Heat shock protein beta-6 protein (1-160AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A