

ID4 Antibody, Biotin conjugated



PACO51625

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcriptional regulator (lacking a basic DNA binding domain) which negatively regulates the basic helix-loop-helix (bHLH) transcription factors by forming heterodimers and inhibiting their DNA binding and transcriptional activity. Implicated in regulating a variety of cellular processes, including cellular growth, senescence, differentiation, apoptosis, angiogenesis, and neoplastic transformation.

Gene ID:

ID4

Uniprot

P47928

Synonyms:

DNA-binding protein inhibitor ID-4 (Class B basic helix-loop-helix protein 27) (bHLHb27) (Inhibitor of DNA binding 4) (Inhibitor of differentiation 4), ID4, BHLHB27

Immunogen:

Recombinant Human DNA-binding protein inhibitor ID-4 protein (1-161AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A