

KIF2A Antibody, HRP conjugated



PACO51675

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plus end-directed microtubule-dependent motor required for normal brain development. May regulate microtubule dynamics during axonal growth. Required for normal progression through mitosis. Required for normal congress of chromosomes at the metaphase plate. Required for normal spindle dynamics during mitosis. Promotes spindle turnover. Implicated in formation of bipolar mitotic spindles. Has microtubule depolymerization activity.

Gene ID:

KIF2A

Uniprot

O00139

Synonyms:

Kinesin-like protein KIF2A (Kinesin-2) (hK2), KIF2A, KIF2 KNS2

Immunogen:

Recombinant Human Kinesin-like protein KIF2A protein (5-308AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A