

MGST2 Antibody, HRP conjugated



PACO51707

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Can catalyze the production of LTC4 from LTA4 and reduced glutathione. Can catalyze the conjugation of 1-chloro-2,4-dinitrobenzene with reduced glutathione.

Gene ID:

MGST2

Uniprot

Q99735

Synonyms:

Microsomal glutathione S-transferase 2 (Microsomal GST-2) (EC 2.5.1.18) (Microsomal GST-II), MGST2, GST2

Immunogen:

Recombinant Human Microsomal glutathione S-transferase 2 protein (27-57AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A