

MYOG Antibody, Biotin conjugated



PACO51733

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

MYOG

Uniprot

P15173

Synonyms:

Myogenin (Class C basic helix-loop-helix protein 3) (bHLHc3) (Myogenic factor 4) (Myf-4), MYOG, BHLHC3 MYF4

Immunogen:

Recombinant Human Myogenin protein (1-224AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A