

PNLIPRP2 Antibody, HRP conjugated



PACO51823

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Lipase with broad substrate specificity. Can hydrolyze both phospholipids and galactolipids. Acts preferentially on monoglycerides, phospholipids and galactolipids. Contributes to milk fat hydrolysis.

Gene ID:

PNLIPRP2

Uniprot

P54317

Synonyms:

Pancreatic lipase-related protein 2 (PL-RP2) (EC 3.1.1.26) (EC 3.1.1.3) (Galactolipase), PNLIPRP2, PLRP2

Immunogen:

Recombinant Human Pancreatic lipase-related protein 2 protein (305-460AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A