

SCTR Antibody, Biotin conjugated



PACO51941

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

This is a receptor for secretin. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.

Gene ID:

SCTR

Uniprot

P47872

Synonyms:

Secretin receptor (SCT-R), SCTR

Immunogen:

Recombinant Human Secretin receptor protein (51-135AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A