

SLITRK5 Antibody



PACO51966

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200, IF:1:50-1:200

Protein Background:

Suppresses neurite outgrowth.

Gene ID:

SLITRK5

Uniprot

O94991

Synonyms:

SLIT and NTRK-like protein 5 (Leucine-rich repeat-containing protein 11), SLITRK5, KIAA0918 LRRC11


Immunogen:

Recombinant Human SLIT and NTRK-like protein 5 protein (617-890AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4


Product Images



Western Blot. Positive WB detected in: Mouse brain tissue. All lanes: SLITRK5 antibody at 2 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 108, 80 kDa. Observed band size: 108 kDa.



Immunocytochemistry analysis of human lung cancer using PACO51966 at dilution of 1:100.



Immunohistochemistry analysis of human testis tissue using PACO51966 at dilution of 1:100.