

# SRR Antibody, Biotin conjugated



PACO51981

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Catalyzes the synthesis of D-serine from L-serine. D-serine is a key coagonist with glutamate at NMDA receptors. Has dehydratase activity towards both L-serine and D-serine.

**Gene ID:**

SRR

**Uniprot**

Q9GZT4

**Synonyms:**

Serine racemase (EC 5.1.1.18) (D-serine ammonia-lyase) (D-serine dehydratase) (EC 4.3.1.18) (L-serine ammonia-lyase) (L-serine dehydratase) (EC 4.3.1.17), SRR

**Immunogen:**

Recombinant Human Serine racemase protein (1-340AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A