

TMEM11 Antibody



PACO52002

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Plays a role in mitochondrial morphogenesis.

Gene ID:

TMEM11

Uniprot

P17152

Synonyms:

Transmembrane protein 11, mitochondrial (Protein PM1) (Protein PMI), TMEM11, C17orf35 PM1


Immunogen:

Recombinant Human Transmembrane protein 11, mitochondrial protein (1-83AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO52002 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human lymph node tissue using PACO52002 at dilution of 1:100.