

TXK Antibody, Biotin conjugated



PACO52029

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TXK

Uniprot

P42681

Synonyms:

Tyrosine-protein kinase TXK (EC 2.7.10.2) (Protein-tyrosine kinase 4) (Resting lymphocyte kinase), TXK, PTK4 RLK

Immunogen:

Recombinant Human Tyrosine-protein kinase TXK protein (32-83AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A