

PACO52030

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Probably involved in intracellular signal transduction by being involved in the initiation of type I IFN signaling. Phosphorylates the interferon-alpha/beta receptor alpha chain.

Gene ID:

TYK2

Uniprot

P29597

Synonyms:

Non-receptor tyrosine-protein kinase TYK2 (EC 2.7.10.2), TYK2


Immunogen:

Recombinant Human Non-receptor tyrosine-protein kinase TYK2 protein (286-372AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO52030 at dilution of 1:100.



Immunofluorescent analysis of A549 cells using PACO52030 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).