

# TYK2 Antibody, Biotin conjugated



PACO52033

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probably involved in intracellular signal transduction by being involved in the initiation of type I IFN signaling. Phosphorylates the interferon-alpha/beta receptor alpha chain.

**Gene ID:**

TYK2

**Uniprot**

P29597

**Synonyms:**

Non-receptor tyrosine-protein kinase TYK2 (EC 2.7.10.2), TYK2

**Immunogen:**

Recombinant Human Non-receptor tyrosine-protein kinase TYK2 protein (286-372AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A