

ALCAM Antibody, HRP conjugated



PACO52079

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

ALCAM

Uniprot

Q13740

Synonyms:

CD166 antigen (Activated leukocyte cell adhesion molecule) (CD antigen CD166), ALCAM, MEMD

Immunogen:

Recombinant Human CD166 antigen protein (28-516AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A