

GAD2 Antibody, Biotin conjugated



PACO52089

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the production of GABA.

Gene ID:

GAD2

Uniprot

Q05329

Synonyms:

Glutamate decarboxylase 2 (EC 4.1.1.15) (65 kDa glutamic acid, decarboxylase) (GAD-65) (Glutamate decarboxylase 65 kDa isoform), GAD2, GAD65

Immunogen:

Recombinant Human Glutamate decarboxylase 2 protein (1-82AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A