

sox10 Antibody, HRP conjugated



PACO52111

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

enteric nervous system development, glial cell differentiation, inner ear development, melanocyte differentiation, negative regulation of neurogenesis, neural crest cell development, neural crest cell migration, olfactory bulb development, oligodendrocyte differentiation, otic vesicle morphogenesis.

Gene ID:

sox10

Uniprot

Q90XD1

Synonyms:

Transcription factor SOX-10, Sox10

Immunogen:

Recombinant Zebrafish HMG-box transcription factor Sox10 protein (1-485AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A