

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

nucleus, DNA binding, transcription regulatory region DNA binding, angiogenesis, convergent extension involved in nephron morphogenesis, determination of left/right symmetry, embryonic heart tube morphogenesis, erythrocyte development, erythrocyte differentiation, heart looping.

Gene ID:

tbx16

Uniprot

Q504I8

Synonyms:

Immunogen:

Recombinant Zebrafish T-box gene 16 protein (1-470AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A