

PACO52122

---

## Product Information

**Size:**

50ug

**Reactivity:**

Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

anterior/posterior pattern specification, Notch signaling pathway, regulation of neurogenesis, regulation of neuron differentiation, somite specification, thalamus development.

**Gene ID:**

her6

**Uniprot**

Q6PBX3

**Synonyms:**

**Immunogen:**

Recombinant Zebrafish Hairy-related 6 protein (1-270AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A