

Larval cuticle protein A1A Antibody, Biotin conjugated

PACO52133



Product Information

Size:

50ug

Reactivity:

Tenebrio molitor

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Component of the cuticle of the larva of Tenebrio molitor.

Gene ID:**Uniprot**

P80681

Synonyms:

Larval cuticle protein A1A (TM-LCP A1A) (TM-A1A)

Immunogen:

Recombinant Tenebrio molitor Larval cuticle protein A1A protein (1-174AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A