

Beta-lactamase inhibitory protein Antibody



PACO52166

Product Information

Size:

50ug

Reactivity:

Streptomyces clavuligerus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Inhibits a wide variety of beta lactamases.

Gene ID:**Uniprot**

P35804

Synonyms:

Beta-lactamase inhibitory protein (BLIP)


Immunogen:

Recombinant Streptomyces clavuligerus Beta-lactamase inhibitory protein (37-201AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: Beta-lactamase inhibitory protein antibody at 4µg/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 33 kDa. Observed band size: 33 kDa.