

PACO52190

Product Information

Size:

50ug

Reactivity:

Bacillus subtilis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

pbpD

Uniprot

P40750

Synonyms:

Penicillin-binding protein 4 (PBP 4) [Includes: Penicillin-insensitive transglycosylase (EC 2.4.1.129) (Peptidoglycan TGase); Penicillin-sensitive transpeptidase (EC 3.4.16.4) (DD-transpeptidase)], pbpD

Immunogen:

Recombinant Bacillus subtilis Penicillin-binding protein 4 protein (213-450AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A