

**Probable terminase, ATPase subunit Antibody,
HRP conjugated
PACO52211**



Product Information

Size:

50ug

Protein Background:

Reactivity:

Haemophilus phage HP1

Gene ID:

Source:

Rabbit

Uniprot

P51718

Isotype:

IgG

Synonyms:

Probable terminase, ATPase subunit (ORF16)

Applications:

ELISA

Immunogen:

Recombinant Haemophilus phage HP1 Probable terminase, ATPase subunit protein (1-190AA).

Recommended dilutions:

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A