

X Antibody, Biotin conjugated



PACO52233

Product Information

Size:

50ug

Protein Background:**Reactivity:**

Hepatitis B virus genotype D subtype ayw

Gene ID:

X

Source:

Rabbit

Uniprot

P03165

Isotype:

IgG

Synonyms:

Protein X (HBx) (Peptide X) (pX), X

Applications:

ELISA

Immunogen:

Recombinant Hepatitis B virus genotype D subtype ayw Protein X protein (1-154AA).

Recommended dilutions:**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A