

Anionic trypsin Antibody



PACO52234

Product Information

Size:

50ug

Reactivity:

Dog

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

extracellular region, extracellular space, calcium ion binding, serine-type endopeptidase activity, collagen catabolic process, digestion, proteolysis.

Gene ID:

Uniprot

P06872

Synonyms:

Anionic trypsin (EC 3.4.21.4)

Immunogen:

Recombinant Dog Anionic trypsin protein (24-247AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A