

Lysozyme C-1 Antibody



PACO52246

Product Information

Size:

50ug

Reactivity:

Anas platyrhynchos

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Lysozymes have primarily a bacteriolytic function; those in tissues and body fluids are associated with the monocyte-macrophage system and enhance the activity of immunoagents.

Gene ID:**Uniprot**

P00705

Synonyms:

Lysozyme C-1 (EC 3.2.1.17) (1,4-beta-N-acetylmuramidase C)


Immunogen:

Recombinant Anas platyrhynchos Lysozyme C-1 protein (19-147AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: Lysozyme C-1 antibody at 2 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 30 kDa. Observed band size: 30 kDa.