

PACO52262

Product Information

Size:

50ug

Reactivity:

Clostridium botulinum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cleaves peptides in various proteins in a process that requires ATP hydrolysis. Has a chymotrypsin-like activity. Plays a major role in the degradation of misfolded proteins.

Gene ID:

clpP

Uniprot

A5I6W1

Synonyms:

ATP-dependent Clp protease proteolytic subunit (EC 3.4.21.92) (Endopeptidase Clp), clpP

Immunogen:

Recombinant Clostridium botulinum ATP-dependent Clp protease proteolytic subunit protein (1-194AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A