

NEUROD1 Antibody, FITC conjugated



PACO52316

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

NEUROD1

Uniprot

Q13562

Synonyms:

Neurogenic differentiation factor 1 (NeuroD) (NeuroD1) (Class A basic helix-loop-helix protein 3) (bHLHa3), NEUROD1, BHLHA3 NEUROD

Immunogen:

Recombinant Human Neurogenic differentiation factor 1 protein (154-356AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A