

# CEBPE Antibody, Biotin conjugated



PACO52321

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

C/EBP are DNA-binding proteins that recognize two different motifs: the CCAAT homology common to many promoters and the enhanced core homology common to many enhancers.

**Gene ID:**

CEBPE

**Uniprot**

Q15744

**Synonyms:**

CCAAT/enhancer-binding protein epsilon (C/EBP epsilon), CEBPE

**Immunogen:**

Recombinant Human CCAAT/enhancer-binding protein epsilon protein (1-200AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A