

MAPRE2 Antibody, HRP conjugated



PACO52343

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in microtubule polymerization, and spindle function by stabilizing microtubules and anchoring them at centrosomes. May play a role in cell migration.

Gene ID:

MAPRE2

Uniprot

Q15555

Synonyms:

Microtubule-associated protein RP/EB family member 2 (APC-binding protein EB2) (End-binding protein 2) (EB2), MAPRE2, RP1

Immunogen:

Recombinant Human Microtubule-associated protein RP/EB family member 2 protein (159-240AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A