

# MAPRE2 Antibody, Biotin conjugated



PACO52345

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May be involved in microtubule polymerization, and spindle function by stabilizing microtubules and anchoring them at centrosomes. May play a role in cell migration.

**Gene ID:**

MAPRE2

**Uniprot**

Q15555

**Synonyms:**

Microtubule-associated protein RP/EB family member 2 (APC-binding protein EB2) (End-binding protein 2) (EB2), MAPRE2, RP1

**Immunogen:**

Recombinant Human Microtubule-associated protein RP/EB family member 2 protein (159-240AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A