

TNFAIP1 Antibody, Biotin conjugated



PACO52389

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Substrate-specific adapter of a BCR (BTB-CUL3-RBX1) E3 ubiquitin-protein ligase complex involved in regulation of cytoskeleton structure. The BCR(BACURD2) E3 ubiquitin ligase complex mediates the ubiquitination of RHOA, leading to its degradation by the proteasome, thereby regulating the actin cytoskeleton and cell migration. Its interaction with RHOB may regulate apoptosis. May enhance the PCNA-dependent DNA polymerase delta activity.

Gene ID:

TNFAIP1

Uniprot

Q13829

Synonyms:

BTB/POZ domain-containing adapter for CUL3-mediated RhoA degradation protein 2 (hBACURD2) (BTB/POZ domain-containing protein TNFAIP1) (Protein B12) (Tumor necrosis factor, alpha-induced protein 1, endothelial), TNFAIP1, BACURD2 EDP1

Immunogen:

Recombinant Human BTB/POZ domain-containing adapter for CUL3-mediated RhoA degradation protein 2 protein (233-316AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A