

LRRK1 Antibody, FITC conjugated



PACO52420

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytoplasm, mitochondrion, identical protein binding, positive regulation of canonical Wnt signaling pathway.

Gene ID:

LRRK1

Uniprot

Q38SD2

Synonyms:

Leucine-rich repeat serine/threonine-protein kinase 1 (EC 2.7.11.1), LRRK1, KIAA1790

Immunogen:

Recombinant Human Leucine-rich repeat serine/threonine-protein kinase 1 protein (1-68AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A