

# LRRK1 Antibody, Biotin conjugated



PACO52421

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

cytoplasm, mitochondrion, identical protein binding, positive regulation of canonical Wnt signaling pathway.

**Gene ID:**

LRRK1

**Uniprot**

Q38SD2

**Synonyms:**

Leucine-rich repeat serine/threonine-protein kinase 1 (EC 2.7.11.1), LRRK1, KIAA1790

**Immunogen:**

Recombinant Human Leucine-rich repeat serine/threonine-protein kinase 1 protein (1-68AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A