

# SLC24A5 Antibody, Biotin conjugated



PACO52473

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Cation exchanger involved in pigmentation, possibly by participating in ion transport in melanosomes. Predominant sodium-Calcium exchanger in melanocytes. Probably transports 1 Ca(2+) and 1 K(+) to the melanosome in exchange for 4 cytoplasmic Na(+).

**Gene ID:**

SLC24A5

**Uniprot**

Q71RS6

**Synonyms:**

Sodium/potassium/calcium exchanger 5 (Na(+)/K(+)/Ca(2+)-exchange protein 5) (Solute carrier family 24 member 5), SLC24A5, JSX NCKX5

**Immunogen:**

Recombinant Human Sodium/potassium/calcium exchanger 5 protein (240-298AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A