

NPSR1 Antibody, FITC conjugated



PACO52480

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

G-protein coupled receptor for neuropeptide S (NPS). Promotes mobilization of intracellular Ca(2+) stores. Inhibits cell growth in response to NPS binding. Involved in pathogenesis of asthma and other IgE-mediated diseases.

Gene ID:

NPSR1

Uniprot

Q6W5P4

Synonyms:

Neuropeptide S receptor (G-protein coupled receptor 154) (G-protein coupled receptor PGR14) (G-protein coupled receptor for asthma susceptibility), NPSR1, GPR154 GPRA PGR14

Immunogen:

Recombinant Human Neuropeptide S receptor protein (1-52AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A