

# GPIHBP1 Antibody, HRP conjugated



PACO52523

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Plays a key role in the lipolytic processing of chylomicrons. Required for the transport of lipoprotein lipase LPL into the capillary lumen.

**Gene ID:**

GPIHBP1

**Uniprot**

Q8IV16

**Synonyms:**

Glycosylphosphatidylinositol-anchored high density lipoprotein-binding protein 1 (GPI-HBP1) (GPI-anchored HDL-binding protein 1) (High density lipoprotein-binding protein 1), GPIHBP1, HBP1

**Immunogen:**

Recombinant Human Glycosylphosphatidylinositol-anchored high density lipoprotein-binding protein 1 protein (21-151AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A