

# ELMOD2 Antibody, HRP conjugated



PACO52527

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Acts as a GTPase-activating protein (GAP) toward guanine nucleotide exchange factors like ARL2, ARL3, ARF1 and ARF6, but not for GTPases outside the Arf family. Regulates IFN-related antiviral responses.

**Gene ID:**

ELMOD2

**Uniprot**

Q8IZ81

**Synonyms:**

ELMO domain-containing protein 2, ELMOD2

**Immunogen:**

Recombinant Human ELMO domain-containing protein 2 protein (82-200AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A