

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Probable flavoprotein which may function in endoplasmic reticulum associated degradation (ERAD). May bind non-native proteins in the endoplasmic reticulum and target them to the ubiquitination machinery for subsequent degradation.

**Gene ID:**

FOXRED2

**Uniprot**

Q8IWF2

**Synonyms:**

FAD-dependent oxidoreductase domain-containing protein 2 (Endoplasmic reticulum flavoprotein associated with degradation), FOXRED2, ERFAD

**Immunogen:**


Recombinant Human FAD-dependent oxidoreductase domain-containing protein 2 protein (501-684AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human brain tissue using PACO52534 at dilution of 1:100.