

VPS39 Antibody, FITC conjugated



PACO52580

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Regulator of TGF-beta/activin signaling, inhibiting SMAD3- and activating SMAD2-dependent transcription. Acts by interfering with SMAD3/SMAD4 complex formation, this would lead to inhibition of SMAD3-dependent transcription and relieve SMAD3 inhibition of SMAD2-dependent promoters, thus increasing SMAD2-dependent transcription. Does not affect TGF-beta-induced SMAD2 or SMAD3 phosphorylation, nor SMAD2/SMAD4 complex formation.

Gene ID:

VPS39

Uniprot

Q96JC1

Synonyms:

Vam6/Vps39-like protein (TRAP1-like protein) (hVam6p), VPS39, KIAA0770 TLP VAM6

Immunogen:

Recombinant Human Vam6/Vps39-like protein (263-490AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A