

# GABPB2 Antibody, FITC conjugated



PACO52648

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May function as transcription factor capable of interacting with purine rich repeats (GA repeats).

**Gene ID:**

GABPB2

**Uniprot**

Q8TAK5

**Synonyms:**

GA-binding protein subunit beta-2 (GABP subunit beta-2) (GABPB-2), GABPB2

**Immunogen:**

Recombinant Human GA-binding protein subunit beta-2 protein (163-440AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A