

Loxl4 Antibody



PACO52658

Product Information

Size:

50ug

Reactivity:

Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

May modulate the formation of a collagenous extracellular matrix.

Gene ID:

Loxl4

Uniprot

Q924C6

Synonyms:

Lysyl oxidase homolog 4 (EC 1.4.3) (Lysyl oxidase-like protein 4) (Lysyl oxidase-related protein C), Loxl4, Loxc


Immunogen:

Recombinant Mouse Lysyl oxidase homolog 4 protein (26-757AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: Rat brain tissue. All lanes: Loxl4 antibody at 4 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 85 kDa. Observed band size: 85 kDa.