

# TRIM54 Antibody, Biotin conjugated



PACO52737

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May bind and stabilize microtubules during myotubes formation.

**Gene ID:**

TRIM54

**Uniprot**

Q9BYV2

**Synonyms:**

Tripartite motif-containing protein 54 (Muscle-specific RING finger protein) (MuRF)  
(Muscle-specific RING finger protein 3) (MuRF-3) (MuRF3) (RING finger protein 30),  
TRIM54, MURF MURF3 RNF30

**Immunogen:**

Recombinant Human Tripartite motif-containing protein 54 protein (60-358AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A