

TRIM5 Antibody, HRP conjugated



PACO52799

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TRIM5

Uniprot

Q9C035

Synonyms:

Tripartite motif-containing protein 5 (EC 2.3.2.27) (RING finger protein 88) (RING-type E3 ubiquitin transferase TRIM5), TRIM5, RNF88

Immunogen:

Recombinant Human Tripartite motif-containing protein 5 protein (1-493AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A