

GPAM Antibody, Biotin conjugated



PACO52805

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Esterifies acyl-group from acyl-ACP to the sn-1 position of glycerol-3-phosphate, an essential step in glycerolipid biosynthesis.

Gene ID:

GPAM

Uniprot

Q9HCL2

Synonyms:

Glycerol-3-phosphate acyltransferase 1, mitochondrial (GPAT-1) (EC 2.3.1.15), GPAM, GPAT1 KIAA1560

Immunogen:

Recombinant Human Glycerol-3-phosphate acyltransferase 1, mitochondrial protein (281-459AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A