

FARSB Antibody, HRP conjugated



PACO52827

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytoplasm, cytosol, membrane, phenylalanine-tRNA ligase complex, phenylalanine-tRNA ligase activity, phenylalanyl-tRNA aminoacylation, protein heterotetramerization, translation, tRNA aminoacylation for protein translation.

Gene ID:

FARSB

Uniprot

Q9NSD9

Synonyms:

Phenylalanine--tRNA ligase beta subunit (EC 6.1.1.20) (Phenylalanyl-tRNA synthetase beta subunit) (PheRS), FARSB, FARSLB FRSB

Immunogen:

Recombinant Human Phenylalanine--tRNA ligase β subunit protein (1-299AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A