

WNT10A Antibody



PACO52846

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule important in CNS development. Is likely to signal over only few cell diameters.

Gene ID:

WNT10A

Uniprot

Q9GZT5

Synonyms:

Protein Wnt-10a, WNT10A


Immunogen:

Recombinant Human Protein Wnt-10a protein (160-372AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO52846 at dilution of 1:100.