

WNT10A Antibody, FITC conjugated



PACO52848

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule important in CNS development. Is likely to signal over only few cell diameters.

Gene ID:

WNT10A

Uniprot

Q9GZT5

Synonyms:

Protein Wnt-10a, WNT10A

Immunogen:

Recombinant Human Protein Wnt-10a protein (160-372AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A