

mib1 Antibody



PACO52914

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

mib1

Uniprot

H9GYU5

Synonyms:

E3 ubiquitin-protein ligase mib1 (EC 2.3.2.27) (Protein mind bomb) (RING-type E3 ubiquitin transferase mib1), mib1, mib

Immunogen:

Recombinant Zebrafish E3 ubiquitin-protein ligase mib1 protein (1-274AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A