

DSG3 Antibody



PACO52946

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Component of intercellular desmosome junctions. Involved in the interaction of plaque proteins and intermediate filaments mediating cell-cell adhesion.

Gene ID:

DSG3

Uniprot

P32926

Synonyms:

Desmoglein-3 (130 kDa pemphigus vulgaris antigen) (PVA) (Cadherin family member 6), DSG3, CDHF6


Immunogen:

Recombinant Human Desmoglein-3 protein (50-615AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: Mouse brain tissue. All lanes: DSG3 antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. predicted band size: 108 kDa. observed band size: 108 kDa.